


Feature												
Model	CMS50-L	MD300CP	CMS50-DL	CMS50-D	MD300-C1	Dive 18705	MD300-C23	MD300-C53	MD300-C634	CMS50-D+	CMS50-E	M-100
Style	Fingertip	Fingertip	Fingertip	Fingertip	Fingertip	Fingertip	Fingertip	Pediatric Finger	Professional Finger	Fingertip	Fingertip	Fingertip
Display	Red LED	Red LED	Red LED	Color OLED	Red LED	Color OLED	Color OLED	Color OLED	Color OLED	Color OLED	Color OLED	Color TFT
Directional Display	1 way	1 way	1 way	6 way	1 way	2 way	6 way	6 way	6 way	6 way	4 way	1 way
Pulse Bar	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Plethysmograph Wave				<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Memory										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Pulse Tone										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
SPO2 Alarm										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
PR Alarm										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
SPO2 Software										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Spot Checking	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Short Term Monitoring				<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Long Term Monitoring							<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Sleep Studies										<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Lanyard Cord	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Carry Case Included								<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Soft Rubber Boot		OPTIONAL	OPTIONAL	OPTIONAL	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
Power	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	2 'AAA'	Lithium / AC / USB	1 'AA'
Battery Life	24 hrs / 2000 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 2500 spot checks	30 hrs / 12 hrs (rec.)	18 hrs continuous	30 hrs Continuous
Wireless Option											<input checked="" type="checkbox"/>	
Warranty	90 Day	1 Year	1 Year	1 Year	1 Year	1 Year	2 Years	2 Years	2 Years	1 Year	1 Year	5 Years
Quality Grade	Entry Level	Good	Good	Good	Very Good	Very Good	Excellent	Excellent	Professional	Good	Excellent	Hospital
Price	\$37.95	\$38.95	\$39.95	\$59.95	\$64.95	\$69.95	79.95	\$79.95	99.95	\$119.95	\$139.95	\$224.95